## HIV-1/HIV-2 Qualitative RNA

Order Name: HIV-1/2 RNA Test Number: 5194321 Revision Date: 12/12/2022

TEST NAME		METHO	DDOLOGY	LOINC CODE
HIV-1/HIV-2 Qualitative RNA		Polym	Polymerase Chain Reaction/Real Time	
SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (650 uL)	Plasma	EDTA (Lavender Top)	Frozen
Alternate 1	1 mL (650 uL)	Serum	Clot Activator SST	Frozen
Instructions	<ul> <li>Specimen Type: Lavender-Top (EDTA) Tube, BD SST Tube, BD Vacutainer PPT Plasma Preparation Tube</li> <li>Specimen Storage: FROZEN Sample may be stored for 24 hours at room temperature. It is recommended that serum or plasma be refrigerated or frozen for shipment.</li> <li>Specimen Collection: Collect specimen in lavender-top (EDTA), gel-barrier (SST tube) or a PPT Plasma Preparation Tube. Whole blood specimens stored at room temperature should be centrifuged within 24 hours of collection. Transfer serum or plasma to a polypropylene screw-cap tube for shipment.</li> <li>Specimen Stability: Ambient: 24 hours; Refrigerated : 7 days; Frozen: 6 weeks (stability provided by manufacturer or literature reference.)</li> <li>Causes for Rejection: Specimen grossly hemolyzed; specimen received at room temperature.</li> </ul>			

GENERAL INFORMATION		
1 - 2 days		
139825		
87538, 87535		
Reference Lab		