## Labcorp Oklahoma, Inc. Test Directory

## **Coronary Artery Disease Panel**

Order Name: Coronary Art DP

Test Number: 5194955 Revision Date: 03/21/2023

TEST NAME		METI	HODOLOGY	LOINC CODE	
Coronary Artery Disease Panel		Polyi	Polymerase Chain Reaction		
SPECIMEN REQUIREM	MENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment	
Preferred	8.5 mL (3 mL)	Whole Blood	ACD Solution A or B (Yellow Top)	Room Temperature	
Alternate 1	8.5 mL (3 mL)	Whole Blood	EDTA (Lavender Top)	Room Temperature	
Alternate 2	1	Saliva	Oragene Dx saliva kit	Room Temperature	
Alternate 3	1	Buccal swab	PurFlock buccal swab kit	Room Temperature	
Instructions	Specimen Type: Whole blood or PurFlock buccal swab kit or Oragene Dx saliva kit  Specimen Volume: 8.5 mL whole blood or PurFlock buccal swab kit or Oragene Dx saliva kit  Mininum Volume: 3 mL whole blood or PurFlock buccal swab kit or Oragene Dx saliva kit  Collection: Standard phlebotomy. Follow PurFlock buccal swab kit or Oragene Dx 500 saliva kit collection instructions. Do not eat, drink, smoke, or chew gum 30 min prior to collection.  Specimen Storage: Maintain specimen at room temperature or refrigerate at 4C Do not freeze.  ?Special Instructions: In cases in which there is a known variant documented in the family, the physician may prefer to order Targeted Variant				

GENERAL INFORMATION			
Expected TAT	14 - 21 days In some cases, additional time may be required for confirmatory or reflex tests.		
Clinical Use	This test includes the following genes: ABCA1, ABCG5, ABCG8, APOA1, APOA5, APOB (556 bp of exon 26), APOC3, LDLR, LDLRAP1, LPL, PCSK9 and PON2.		
Performing Labcorp Test Code	482243		
Notes	Labcorp Test Code: 482243		
Lab Section	Reference Lab		

Analysis, test code 5194970. Test orders must include an attestation that the provider has the patient's informed consent for genetic testing.

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2025