

Proinsulin

Order Name: **PROINSULIN**

Test Number: 3655950

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Proinsulin	Enzyme immunoassay (EIA)	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.4)	Serum	Clot Activator SST	Frozen
Alternate 1	1 mL (0.4)	Serum	Clot Activator (Red Top, No-Gel)	Frozen
Instructions	<p><b>Notes:</b> 0.4 mL (Note: This volume Does NOT allow for repeat testing.)</p> <p><b>Specimen Type:</b> Red-top tube or gel-barrier tube</p> <p><b>Specimen Storage:</b> <b>FREEZE</b></p> <p><b>Specimen Collection:</b> Transfer the serum into a LabCorp PP transpak purple tube with screw cap (LabCorp ID 49482). Freeze immediately and maintain frozen until tested. To avoid delays in turnaround time when requesting multiple test on frozen samples, PLEASE SUBMIT SEPARATE FROZEN SPECIMENS FOR EACH TEST REQUESTED.</p> <p><b>Specimen Stability:</b> Ambient: Unstable, Refrigerated : Unstable, Frozen: 12 days</p>			

GENERAL INFORMATION	
Expected TAT	3 - 5 days
Performing Labcorp Test Code	140533
Notes	Labcorp Test Code: 140533
CPT Code(s)	84206
Lab Section	Reference Lab