Labcorp Oklahoma, Inc. Test Directory

Felbamate (Felbatol)

Order Name: **FELBAMATE**Test Number: 3657650
Revision Date: 12/12/2022

TEST NAME			METHODOLOGY	LOINC CODE
Felbamate (Felbato	ol)		Liquid Chromatography/Tandem Mass Spectrom	netry 6899-9
SPECIMEN REQUI	REMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (0.6 mL)	Serum	Clot Activator (Red Top, No-Gel)	Room Temperature
Alternate 1	2 mL (0.6 mL)	Plasma	EDTA (Lavender Top)	Room Temperature
Alternate 2	2 mL (0.6 mL)	Plasma	Sodium Heparin (Green Top / No- Gel)	Room Temperature
Instructions	Specimen Type: Red-top tube, lavender-top (EDTA) tube OR green-top (heparin) tube. DO NOT USE A GEL-BARRIER TUBE. The use of gelbarrier tubes is not recommended due to slow absorption of the drug by the gel. Depending on the specimen volume and storage time, the decrease in drug level due to absorption may be clinically significant. Specimen Storage: Room Temperature Specimen Collection: Transfer separated serum or plasma to a plastic transport tube. Specimen Stability: Ambient: 14 days, Refrigerated: 14 days, Frozen: 14 days			

GENERAL INFORMATION			
Expected TAT	2 - 4 days one working day.		
Clinical Use	Felbamate is an antiepileptic drug used to treat patients with a variety of types of seizures. Therapeutic drug monitoring is useful to optimize dose and avoid toxicity		
Performing Labcorp Test Code	716530		
Notes	Labcorp Test Code: 716530		
CPT Code(s)	80167		
Lab Section	Reference Lab		

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2026