## Labcorp Oklahoma, Inc. Test Directory

## HIV-1/HIV-2 Qualitative RNA

Order Name: HIV-1/2 RNA
Test Number: 5194321
Revision Date: 12/12/2022

TEST NAME			METHODOLOGY	LOINC CODE
HIV-1/HIV-2 Qualitative RNA			Polymerase Chain Reaction/Real Time	
SPECIMEN REQU	JIREMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (650 uL)	Plasma	EDTA (Lavender Top)	Frozen
Alternate 1	1 mL (650 uL)	Serum	Clot Activator SST	Frozen
Instructions	Specimen Type: Lavender-To	Specimen Type: Lavender-Top (EDTA) Tube, BD SST Tube, BD Vacutainer PPT Plasma Preparation Tube		

Specimen Storage: FROZEN Sample may be stored for 24 hours at room temperature. It is recommended that serum or plasma be refrigerated or frozen for shipment.

Specimen Collection: Collect specimen in lavender-top (EDTA), gel-barrier (SST tube) or a PPT Plasma Preparation Tube. Whole blood specimens

specimen Collection: Collect specimen in lavender-top (EDTA), gel-barrier (SST tube) or a PPT Plasma Preparation Tube. Whole blood specimens stored at room temperature should be centrifuged within 24 hours of collection. Transfer serum or plasma to a polypropylene screw-cap tube for shipment.

**Specimen Stability:** Ambient: 24 hours; Refrigerated: 7 days; Frozen: 6 weeks (stability provided by manufacturer or literature reference.) **Causes for Rejection**: Specimen grossly hemolyzed; specimen received at room temperature.

GENERAL INFORMATION		
Expected TAT	1 - 2 days	
Performing Labcorp Test Code	139825	
CPT Code(s)	87538, 87535	
Lab Section	Reference Lab	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2026