

Rotavirus Screen

Order Name: **C ROTA SC**

Test Number: 6100250

Revision Date: 02/02/2016

TEST NAME	METHODOLOGY	LOINC CODE
Rotavirus Screen	Enzyme Immunoassay	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	3mL	Fecal/Stool	Sterile Orange Screwtop Container	Refrigerated
Alternate 1	3mL	Fecal/Stool	Stool Specimen Container, Small for Random Collection (Gray)	Refrigerated
Instructions	Specimen is good for 72 hours refrigerated.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1 Day
Clinical Use	Detects rotavirus in stool
CPT Code(s)	87425
Lab Section	Microbiology