

Bakers Yeast Allergen IgG

Order Name: **YEAST IgG**
Test Number: 3666150
Revision Date: 12/14/2022

| TEST NAME | METHODOLOGY | LOINC CODE |
|---------------------------|-------------|------------|
| Bakers Yeast Allergen IgG | ImmunoCAP | 35538-8 |

SPECIMEN REQUIREMENTS

| | | | | |
|--------------|--|---------------|---------------------------|-------------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 0.3 mL (0.1 mL) | Serum | Clot Activator SST | Room Temperature |
| Instructions | <p>Specimen Type: Red-top tube or gel-barrier tube Specimen Storage: Room temperature Specimen Collection: Not Available Special Instructions: Values obtained with different assay methods should not be used interchangeably in serial testing. It is recommended that only one assay method be used consistently to monitor each patient's course of therapy. Specimen Stability: Ambient: 14 days, Refrigerated : 14 days, Frozen: 3 months</p> | | | |

GENERAL INFORMATION

| | |
|------------------------------|---------------------------|
| Expected TAT | 2 - 4 days |
| Performing Labcorp Test Code | 605066 |
| Notes | Labcorp Test Code: 605066 |
| CPT Code(s) | 86001 |
| Lab Section | Reference Lab |