

Donath-Landsteiner Test

Order Name: Donath-Land

Test Number: 7306200

Revision Date: 09/24/2019

| TEST NAME | METHODOLOGY | LOINC CODE |
|-------------------------|-------------|------------|
| Donath-Landsteiner Test | Hemolysis | |

| SPECIMEN REQUIREMENTS | | | | |
|-----------------------|--|---------------|----------------------------------|-----------------------|
| Specimen | Specimen Volume (min) | Specimen Type | Specimen Container | Transport Environment |
| Preferred | 10 mL (7 mL) | Serum | Clot Activator (Red Top, No-Gel) | Refrigerated |
| Instructions | Blood should be allowed to clot at 37°C, and the serum separated from the clot at this temperature. Once the serum is separated from the cells, the serum can be stored and transported refrigerated. | | | |

| GENERAL INFORMATION | |
|---------------------|--|
| Testing Schedule | Mon-Fri (7:00am - 3:30pm) |
| Expected TAT | 1 Day |
| Clinical Use | Used in the diagnosis of Paroxysmal Cold Hemoglobinuria (PCH). |
| CPT Code(s) | 86940 86941 |
| Lab Section | Blood Bank |