

Zucchini Allergen IgE

Order Name: **Zucchini**

Test Number: **5613527**

Revision Date: **08/21/2017**

TEST NAME	METHODOLOGY	LOINC CODE
Zucchini Allergen IgE	ImmunoCAP	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	0.25 mL (0.1 mL)	Serum	Clot Activator SST	Refrigerated
Instructions	Specimen: Separate serum from cells ASAP or within 2 hours of collection. **Transfer 0.25mL(0.1mL) serum per allergen ordered to a Standard Transport Tube. **Unacceptable Conditions: Hemolyzed, icteric, or lipemic specimens. **Stability After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year			

GENERAL INFORMATION	
Testing Schedule	Sun-Sat
Expected TAT	1-2 days
CPT Code(s)	86003
Lab Section	Reference Lab