## Flow Cytometry on Bone Marrow

Order Name: FLOW BM
Test Number: 5582550
Revision Date: 07/06/2022

| TEST NAME                     |  |               | METHODOLOGY         | LOINC CODE            |  |
|-------------------------------|--|---------------|---------------------|-----------------------|--|
| Flow Cytometry on Bone Marrow |  |               | Flow cytometry      |                       |  |
| SPECIMEN REQU                 | UREMENTS   |               |                     |                       |  |
| Specimen                      | Specimen Volume (min)  | Specimen Type | Specimen Container  | Transport Environment |  |
| Preferred                     | 5 mL (3 mL)  | Bone Marrow   | EDTA (Lavender Top) | Room Temperature      |  |
| Instructions                  | Please deliver to performing laboratory (flow cytometry) ASAP. |               |                     |                       |  |
|                               | Specimen stability: 48hrs refri                                | gerated       |                     |                       |  |

| GENERAL INFORMATION |   |  |
|---------------------|---|--|
| Testing Schedule    | Tue-Sat   |  |
| Expected TAT        | 2 Days  |  |
| Clinical Use        | Assist in the diagnosis of a lymphoproliferative disease. |  |
| CPT Code(s)         | test dependant  |  |
| Lab Section         | Flow Cytometry  |  |
|                     |   |  |

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2026