

Opiates, Screen Urine

Order Name: **OPIATES U**

Test Number: 4301300

Revision Date: 05/28/2021

TEST NAME		METHODOLOGY	LOINC CODE
Opiates, Screen Urine		CEDIA	3879-4

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	30 mL (10 mL)	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Keep refrigerated. Specimen stability: Refrigerated 5 days, Frozen >5 days			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1-2 days
CPT Code(s)	80307
Lab Section	Chemistry