

Wet Prep for Yeast and Trichomonas

Order Name: **Wet Prep**
Test Number: 2915475
Revision Date: 11/15/2021

TEST NAME	METHODOLOGY	LOINC CODE
Vaginal Yeast Examination	Microscopy	32765-0
Trichomonas Microscopic Examination	Microscopy	32766-8

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	See Instructions	See Instructions	Wet prep saline	Room Temperature
Alternate 1	Swab	Swab	Copan eSwab - White (Regular size)	Room Temperature
Instructions	Collect discharge and place in sterile tube with saline or Copan eSwab - White. Deliver specimens to the laboratory for testing as soon as possible, ideally within 2 hours of collection. Trichomonas viability decreases in the specimen overtime. The Trichomonas antigen (2915460) testing is strongly recommended. Specimens arriving 24 hours after collection are unacceptable.			

GENERAL INFORMATION	
Testing Schedule	Daily
Expected TAT	1 Day
Clinical Use	Determines yeast vaginitis and/or trichomoniasis
Notes	This used to be known as: Order Name C WET PR Test Number 6000650
CPT Code(s)	87210
Lab Section	Microbiology