

## Mirtazapine (Remeron), Serum or Plasma

Order Name: **Mirtazapine**

Test Number: 6905779

Revision Date: 10/20/2022

TEST NAME	METHODOLOGY	LOINC CODE
Mirtazapine (Remeron), Serum or Plasma	High-pressure liquid chromatography (HPLC) with fluorometric detection	17283-3

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	2 mL (0.5 mL)	Serum	Clot Activator (Red Top, No-Gel)	Room Temperature

**Instructions**

**Specimen Type:** Red-top tube OR green-top (heparin) tube. Gel-barrier tubes are NOT recommended.

**Specimen Storage:** SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen.

**Specimen Collection:** Serum or plasma should be separated from cells within two hours of venipuncture. Submit serum or plasma in a plastic transport tube.

**Special Instructions:** Testing referred to MedTox Laboratories Inc MX#434

### GENERAL INFORMATION

Expected TAT	10-14 days
CPT Code(s)	80299
Lab Section	Reference Lab