

ToxAssure Select 13 EDI-DIS, Urine

Order Name: **TA 13 EDI-DIS**
Test Number: 5195570
Revision Date: 01/09/2024

TEST NAME	METHODOLOGY	LOINC CODE
ToxAssure Select 13 EDI-DIS, Urine	Liquid Chromatography/Tandem Mass Spectrometry	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	30 mL (15 mL)	Urine, Random	Sterile Urine Container Preservative Free	Room Temperature
Instructions	Specimen: 30mL(15mL) Random Urine Preservative free Urine Transport tube Room Temperature or Refrigerated. Cause for Rejection: Insufficient volume, no ID on container, urine from preservative tube. SUBMISSION/TRANSPORT (<3 days): Room Temperature. For storage beyond three days, specimen should be refrigerated or frozen. Stability: Room Temperature 14 days, Refrigerated 14 days, Frozen 6 weeks.			

GENERAL INFORMATION	
Expected TAT	5-8 days
Performing Labcorp Test Code	700318
CPT Code(s)	80353, 80348, 80361, 80358, 80359, 80347, 80324, 80354, 80365, 80307, 80373, 80372
Lab Section	Reference Lab