

ENT Allergy Panel 2 (Moderate panel)

Order Name: **ENT PNL 2**

Test Number: 6907457

Revision Date: 08/01/2023

TEST NAME	METHODOLOGY	LOINC CODE
ENT Allergy Panel 2 (Moderate panel)	ImmunoCAP	

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	15 mL (7.8 mL)	Serum	Clot Activator (SST or Red No-Gel)	Room Temperature
Instructions	Specimen Type: Red-top tube or gel-barrier tube, Separate Serum from Cells into a screwtop transport container. Stability Requirements: Room temperature 14 days, Refrigerated 14 days, Frozen 3 months. (Freeze/thaw cycles Stable x3)			

GENERAL INFORMATION			
Expected TAT	3 - 5 days		
Notes	Alternaria alternata Allergen	Common Silver Birch Tree Allergen	Mesquite Tree Allergen
	Aspergillus fumigatus Allergen	Cottonwood Tree Allergen	Mountain Juniper Tree Allergen
	Bermuda Grass Allergen	Dandelion Allergen	Mugwort Allergen
	Bipolaris spicifera IgE	Dog Dander Allergen	Mulberry Tree White Allergen
	Box Elder (Maple) Tree Allergen	Dust Mite (D. farinae) Allergen	Pecan, Hickory Tree Allergen
	Cat Dander Allergen	Elm Tree Allergen	Penicillin Chrysogenium Allergen
	Cephalosporium IgE	Epicoccum purpurascens Allergen	Plantain (English), Ribwort Allergen
	Chaetomium globosum IgE	Firebush (Kochia) Allergen	Privet Tree Allergen
	Cladosporium herbarum Allergen	Giant Ragweed Allergen	Ragweed Common Short Allergen
	Cocklebur Allergen	Goosefoot, Lamb's Quarters Allergen	Rough Marshelder Allergen
	Cockroach, American IgE	House Dust Mite Allergen	Saltwort, Russian Thistle Allergen
	Cockroach, German Allergen	Immunoglobulin IgE	Timothy Grass Allergen
	Common Pigweed Allergen	Johnson Grass Allergen	White Pine Tree IgE
	Allergens performed at Labcorp Phoenix		
CPT Code(s)	86003 x39		
Lab Section	Reference Lab		