Eosinophil, Urine

Order Name: **EOS URINE**Test Number: 1001600
Revision Date: 04/29/2019

TEST NAME			METHODOLOGY	LOINC CODE
Eosinophil, Urine			Wright Stain	12210-1
SPECIMEN REQUIR	REMENTS			
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (1 mL)	Urine, Random	Sterile Urine container	Refrigerated
Instructions	Early morning specimens preferred. Keep specimen refrigerated, deliver to lab ASAP. Stability 24 Hours Refrigerated.			

Testing Schedule Daily	
Expected TAT 1 Day	
CPT Code(s) 85999	
Lab Section Hematology	

Service provided by Labcorp Oklahoma, Inc. All Rights Reserved. © 2003 - 2026