

## Alpha-Fetoprotein (AFP), Maternal Tetra Profile

Order Name: **AFP Tetra**  
Test Number: 5194839  
Revision Date: 12/08/2022

TEST NAME	METHODOLOGY	LOINC CODE
Alpha-Fetoprotein (AFP), Maternal Tetra Profile	Chemiluminescence Assay	

### SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (mL)	Specimen Type	Specimen Container	Transport Environment
Preferred	5 mL (3 mL)	Serum	Clot Activator SST	Room Temperature

**Instructions**

**Notes:** 3 mL (Note: This volume Does NOT allow for repeat testing).

**Specimen Type:** Gel-barrier tube

**Specimen Storage:** Room temperature

**Specimen Collection:** Collect in serum separator tube with gel barrier. Allow blood to clot, avoiding hemolysis. Separate serum from cells by centrifugation. Transport spun tube to testing laboratory. Pour off is not advised. Maternal serum specimens must be drawn prior to amniocentesis to avoid contamination with fetal blood.

**Special Instructions:** The following information must be provided: gestational age, date on which the patient was the stated gestational age, how gestational age was determined (LMP, EDD, US), patient's weight, patient's date of birth, patient's race (white, black, other), and insulin-dependent diabetic status. Also indicate relevant patient history (eg, prior neural tube defects, Down syndrome, ultrasound anomalies, or previous maternal serum screening specimen during this pregnancy). Complete information is necessary to interpret the test. Patient information may be provided to the laboratory using The Maternal Prenatal Screening requisition form 0900. Specimens must be collected before amniocentesis.

**Specimen Stability:** Ambient: 7 days, Refrigerated : 14 days, Frozen: 14 days

**Labcorp Oklahoma, Inc.**  
**Test Directory**

**GENERAL INFORMATION**

<b>Expected TAT</b>	2 - 5 days
<b>Clinical Use</b>	Down syndrome screening is offered for gestational ages 15.0 to 21.9 weeks. Open spina bifida screening is offered for gestational ages 15.0 to 23.9 weeks. The optimal gestational age for open spina bifida screening is 16.0 to 18.9 weeks.
<b>Performing Labcorp Test Code</b>	017319
<b>Notes</b>	Labcorp Test Code: 017319

**Prompt Information - (Please provide as many as possible for the best interpretation)**

REQUIRED	PROMPT CODE	MNEMONIC	LONG NAME	RESULT TYPE
Required Prompt	5194810	LC PAWLBS	Weight (lbs)	5 Dig Numeric
	5195221	LC PAWOZS	Weight (Ounces)	2 Dig Numeric
	5194811	LC INSDEP	Insulin Dependant	Y or N
Required Prompt	5194684	LC GESAWK	Gest Age Weeks (##)	2 Dig Numeric
Required Prompt	5194685	LC GESADY	Gest Age Days (#)	1 Numeric (0-6)
	5195222	LC GESADC	Gest Age In Decimal Form	1 Numeric (0-6)
	5194686	LC GESADT	Gest Age Date of Calc (YYYYMMDD)	YYYYMMDD
	5194687	LC GACMTH	Gest Age Calc Method	LMP, US, EDD/EDC
Requested if by LMP	5194688	LC LMDATE	LMP Date (YYYYMMDD)	YYYYMMDD
Requested if by EDD/EDC	5194689	LC EDDATE	EDD/EDC Date (YYYYMMDD)	YYYYMMDD
	5194812	LC NFETUS	Number of Fetuses	1 Dig Numeric (1-9)
Requested if Present	5194690	LC OTHIND	Other Indications (Y/N)	Y or N
	5194691	LC ADINFO	Additional Info	0-20 characters
	5194692	LC PRELEV	Prev Elevated AFP (Y/N)	Y or N
	5194693	LC DONEGG	Donor Egg (Y/N)	Y or N
	5194694	LC EGGDAG	Age of Egg Donor	2 Dig Numeric
	5194695	LC EGGTYP	Type of Egg Donor	S-SELF, N-NON-SELF
	5194696	LC PDONTD	Prior DS/ONTD Scr Current Preg (Y/N)	Y or N
	5194697	LC PFTTST	Prior 1st Trimester Testing (Y/N)	Y or N
	5194698	LC PSTTST	Prior 2nd Trimester Testing (Y/N)	Y or N
	5194813	LC FHONTD	Family Hx NTD	Y or N
	5194814	LC PPRGWD	Prior Pregnancy with DS	Y or N

<b>CPT Code(s)</b>	86336, 84702, 82105, 82677
<b>Lab Section</b>	Reference Lab

