

Metanephrines, Fractionated, Random Urine

Order Name: **METS PHEO UR**

Test Number: 5194542

Revision Date: 12/12/2022

TEST NAME	METHODOLOGY	LOINC CODE
Metanephrines, Fractionated, Random Urine	Liquid Chromatography/Tandem Mass Spectrometry	

SPECIMEN REQUIREMENTS

Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	10 mL (2.5 mL)	Urine, Random	Sterile Screwtop Container	Room Temperature
Instructions	<p>Notes: 2.5 mL (Note: This volume Does NOT allow for repeat testing.) Specimen Type: Random urine transport tube with no preservative or 6N HCl. Specimen Storage: Room Temperature Specimen Collection: Collect urine and transfer to a random urine transport tube. Specimen Stability: Ambient: 7 days, Refrigerated : 14 days, Frozen: 14 days</p>			

GENERAL INFORMATION

Expected TAT	3 - 7 days
Clinical Use	Metanephrines, Pheochromocytoma Evaluation - Diagnose the presence of a pheochromocytoma, neuroblastoma, and ganglioneuroblastoma Methodology : Metanephrines: liquid chromatography/tandem mass spectrometry (LC/MS-MS); creatinine: alkaline picric acid
Performing Labcorp Test Code	004044
Notes	Labcorp Test Code: 004044
CPT Code(s)	83835, 82570
Lab Section	Reference Lab