

Pneumococcal Ab Pnl PPSV23, 23 Serotypes

Order Name: **Pneumo 23**

Test Number: 5197338

Revision Date: 03/16/2026

TEST NAME	METHODOLOGY	LOINC CODE
Pneumococcal Ab Pnl PPSV23, 23 Serotypes	Multiplex fluorescent immunoassay	See Below

SPECIMEN REQUIREMENTS				
Specimen	Specimen Volume (min)	Specimen Type	Specimen Container	Transport Environment
Preferred	1 mL (0.5 mL)	Serum	Clot Activator (SST or Red No-Gel)	Refrigerated
Instructions	Specimen: 1mL(0.5mL) Serum from SST clot tube Blood should be allowed to clot 30-60 minutes, then poured off into a plastic transport tube Storage Instructions: Refrigerated Stability Requirements: Ambient 6 weeks, Refrigerated 6 weeks, Frozen 6 weeks Causes for Rejection: Insufficient sample volume.			

Labcorp Oklahoma, Inc.
Test Directory

GENERAL INFORMATION						
Expected TAT		7-10 Days				
Clinical Use		Pneumococcal antibody panel provides quantitative measurement of serotype-specific IgG antibodies to assess immune response to pneumococcal vaccination, particularly in immunocompromised patients.				
Performing Labcorp Test Code		834968				
Notes						
		RESULT CODE	RESULT NAME	RESULT LONG NAME	LOINC CODE	
		5575653	Pneumo 1 IgG	Pneumococcal Serotype 1 IgG	85955-3	
		5575657	Pneumo 2 IgG	Pneumococcal Serotype 2 IgG	86038-7	
		5575659	Pneumo 3 IgG	Pneumococcal Serotype 3 IgG	86081-7	
		5575661	Pneumo 4 IgG	Pneumococcal Serotype 4* IgG	86108-8	
		5575663	Pneumo 5 IgG	Pneumococcal Serotype 5 IgG	86129-4	
		5575689	Pneumo 6B IgG	Pneumococcal Serotype 6B* IgG	40905-2	
		5575697	Pneumo 7F IgG	Pneumococcal Serotype 7F IgG	40911-0	
		5575667	Pneumo 8 IgG	Pneumococcal Serotype 8 IgG	86148-4	
		5575669	Pneumo 9N IgG	Pneumococcal Serotype 9N IgG	86166-6	
		5575707	Pneumo 9V IgG	Pneumococcal Serotype 9V*, IgG	40926-8	
		5575691	Pneumo 10A IgG	Pneumococcal Serotype 10A IgG	86099-9	
		5575693	Pneumo 11A IgG	Pneumococcal Serotype 11A IgG	86123-7	
		5575671	Pneumo 12F IgG	Pneumococcal Serotype 12F IgG	85974-4	
		5575673	Pneumo 14 IgG	Pneumococcal Serotype 14* IgG	85992-6	
		5575699	Pneumo 15B IgG	Pneumococcal Serotype 15B IgG	47356-1	
		5575677	Pneumo 17F IgG	Pneumococcal Serotype 17F IgG	86010-6	
		5575701	Pneumo 18C IgG	Pneumococcal Serotype 18C* IgG	40913-6	
		5575703	Pneumo 19A IgG	Pneumococcal Serotype 19A IgG	40915-1	
		5575679	Pneumo 19F IgG	Pneumococcal Serotype 19F* IgG	86021-3	
		5197344	Pneumo 20A IgG	Pneumococcal Serotype 20A IgG	n/a	
		5575683	Pneumo 22F IgG	Pneumococcal Serotype 22F IgG	86053-6	
		5575687	Pneumo 23F IgG	Pneumococcal Serotype 23F* IgG	86061-9	
		5575709	Pneumo 33F IgG	Pneumococcal Serotype 33F IgG	47349-6	
		CPT Code(s)		86581		
Lab Section		Reference Lab				